

# THE SCD IMPLANT BUNDLE

Our best bundle. Without any compromise.

Our SCD Implant Bundle now offers quality components, greater compatibility with popular implant systems and the assurance of an all-inclusive price. The new Titanium Nitride (TiN) coated screw ensures reliable retention, better torqueing, and even better results for your patients.

ALL-INCLUSIVE bundle at only **€305**





Screw- or cement-retained PFM crown*	+	Custom milled abutment	+	<b>NEW</b> Titanium Nitride (TiN) coated screw	+	Implant analogue	+	Soft tissue model
--------------------------------------	---	------------------------	---	---	---	------------------	---	-------------------

\*Semi or high precious alloy extra ~Standard Titanium screw available. ^T&Cs apply, see website for details.

Now available with Titanium Nitride TiN coating

NEW



TITANIUM NITRIDE (TiN) COATING <sup>NEW</sup>

Choose the screw that's right for you.  
Titanium | TiN coated

Every implant case is unique. Therefore we are giving you the choice. Select between Titanium or our new Titanium Nitride (TiN) coated screws. The TiN coating is an important upgrade to our product as this special surface treatment improves the performance of the screw.

- The improved, harder TiN surface provides a more stable and reliable connection with the implant.

## Screw- or Cement-Retained PFM Crown-

All our products are manufactured following globally recognised processes and our restorations are made under a quality management system that is certified under ISO 9001:2008 and ISO 13485:2003.



## Custom Milled Abutment

For our custom, titanium milled abutments we use Grade 4 or Grade 5 quality. All components are tested under ISO 14801:2007 and we use CAD/CAM precision milling to ensure high accuracy.



# IMPLANT\* COMPONENT SUMMARY

Brand	System	Platform	Digitek <sup>†</sup> Abutment	Digitek <sup>^</sup> Bar/Bridge	Original Custom Milling Abutment	Torque Ncm	Bridge Angle <sup>‡</sup>	
3i Biomet	Certain (BC)	3.4	✓	✓	—	20	40°	
		4.1	✓	✓	—	20	40°	
		5.0	✓	✓	—	20	40°	
	External Hex (BO)	3.4	✓	✓	—	20	30°	
		4.1	✓	✓	—	20	40°	
		5.0	✓	✓	—	20	40°	
Zimmer	Tapered Screw-Vent (ZT)	3.5	✓	✓	—	30	60°	
		4.5	✓	✓	—	30	60°	
		5.7	✓	✓	—	30	60°	
Dentsply	AstraTech OsseoSpeed (DA)	3.0	✓	✓	—	15	20°	
		3.5/4.0	✓	✓	—	20	20°	
		4.5/5.0	✓	✓	—	25	20°	
	UniAbutment (DAU)	20°	—	✓	—	15	40°	
		45°	—	✓	—	15	60°	
	AstraTech EV (DV)	3.6	✓*	✓	—	25	20°	
		4.2	✓*	✓	—	25	20°	
		4.8	✓*	✓	—	25	20°	
	Frialit/XIVE (DX)	UniAbutment EV (DVU)	Straight	—	✓	—	15	40°
			3.4	✓	✓	—	24	40°
		MP Abutment (DXM)	3.8	✓	✓	—	24	40°
			4.5	✓	✓	—	24	40°
5.5			✓	✓	—	24	40°	
3.8			—	✓	—	24	40°	
Ankylos (DY)		4.5	—	✓	—	24	40°	
		5.5	—	✓	—	24	40°	
		C/X	—	✓	—	15	10°	
		—	✓	✓	—	35	40°	
Straumann	Tissue Level (ST)	NNC	✓	✓	—	35	30°	
		RN	✓	✓	✓*	35	40°	
	Bone Level (SB)	WN	✓	✓	✓*	35	40°	
		1.5 Abutment (STU)	RN	—	✓	—	15	40°
		NC	✓	✓	✓*	35	30°	
Nobel Biocare	Replace (NR)	RC	✓	✓	✓*	35	30°	
		NP	✓	✓	—	35	40°	
		RP	✓	✓	—	35	40°	
		WP	✓	✓	—	35	40°	
		6.0	✓	✓	—	35	40°	
	Active/Conical (NA)	3.0	✓	✓	—	15	30°	
		NP	✓	✓	—	35	40°	
		RP	✓	✓	—	35	40°	
		WP	✓	✓	—	35	40°	
		NP	✓	✓	—	35	30°	
	Branemark (NB)	RP	✓	✓	—	35	40°	
		WP	✓	✓	—	35	40°	
		NP	✓	✓	—	35	40°	
	MUA (NU)	RP	—	✓	—	15	40°	
		WP	—	✓	—	15	40°	
On1 (NO)	NP	—	—	—	35	?		
	RP	—	PLANNED	—	35	?		
	WP	—	—	—	35	?		
Osstem	TS (OT)	Mini	✓	✓	—	20	20°	
		Standard	✓	✓	—	30	20°	

\*Only titanium available

†One-piece design, precision machined in titanium, CoCr

Torque - recommended screw torque

‡Angle, for bridge, max angle between two implants axis

For a complete list of supported implant systems, please visit our website.

# IMPLANT\* COMPONENT SUMMARY

Brand	System	Platform	Digitek <sup>†</sup> Abutment	Digitek <sup>^</sup> Bar/Bridge	Original Custom Milling Abutment	Torque Ncm	Bridge Angle <sup>‡</sup>	
Camlog	Camlog (CC)	3.3	✓*	✓	—	20	40°	
		3.8	✓*	✓	—	20	40°	
		4.3	✓*	✓	—	20	40°	
		5.0	✓*	✓	—	20	40°	
	Conelog (CO)	3.3	✓*	✓	—	20	40°	
		3.8/4.3	✓*	✓	—	20	40°	
Anthogyr	Axiom (AA)	5.0	—	✓	—	20	40°	
		—	✓	✓	—	25	10°	
	Multi Unit (AAU)	Narrow	—	✓	—	15	40°	
Biotech	Kontakt (TK)	Standard	—	✓	—	15	40°	
		—	✓	✓	—	20	10°	
	Conical Abut. (TKU)	Narrow	—	✓	—	20	40°	
AB	AB (AB)	Standard	—	✓	—	20	40°	
		Narrow	✓	✓	—	30	60°	
	Multi Unit (ABU)	P14	—	✓	—	20	40°	
		P16	—	✓	—	20	40°	
BTI	Interna (TI)	4.1	✓	✓	—	35	40°	
		Multi Im (TUI)	Straight	—	✓	—	20	40°
	Angulado	—	✓	—	20	40°		
Medentis	ICX (MEI)	—	✓*	✓	—	30	20°	
		NP	✓*	✓	—	30	40°	
	C1 (MIC)	SP	✓*	✓	—	30	40°	
		WP	✓*	✓	—	30	40°	
		V3	NP	✓*	✓	—	30	40°
	MIS	Multi unit (MIU)	SP	✓*	✓	—	30	40°
			4.8	—	✓	—	25	40°
		Seven (MIS)	NP	✓*	✓	—	30	60°
			SP	✓*	✓	—	30	60°
	Euroteknika	Nat. (EUN)	WP	✓*	✓	—	30	60°
Standard			✓	✓	—	25	20°	
Conical Abut. (EUU)		Tetra	—	✓	—	20	40°	
GlobalD	In-Kone (GLI)	Plural	—	✓	—	20	40°	
		—	✓	✓	—	15	8°	
	Conical Abut. (GLU)	Straight	—	✓	—	15	40°	
Bego	Semados (BES)	Angled	—	✓	—	15	40°	
		3.25/3.75	✓	✓	—	30	60°	
		4.1	✓	✓	—	30	60°	
		4.5	✓	✓	—	30	60°	
Biohorizons	Biohorizons (BHB)	4.5	✓	✓	—	30	60°	
		5.5	✓	✓	—	30	60°	
		3.5	—	—	—	30	60°	
Dentium	Superline/Implantium (DTS)	4.5	—	—	—	30	60°	
		—	—	—	—	30	20°	
	Screw Abutment (DTU)	4.5	—	—	—	25	40°	

\*Only titanium available

†One-piece design, precision machined in titanium, CoCr

Torque - recommended screw torque

‡Angle, for bridge, max angle between two implants axis

For a complete list of supported implant systems, please visit our website.